




ABRASIVES

BROWNIE®/GREENIE®/SUPERGREENIE POLISHERS



- Silicone polishers impregnated with high quality abrasives
- Use **Brownies** for pre-polishing, **Greenies** for polishing & **Supergreenies** for super-polishing
- Use **Brownie/Greenie/Supergreenie Stones for:** amalgams, precious alloys & non-precious alloys

SHAPE CATEGORY	Cup	Mini-Pt.	KN7	PC2	WH6	Floppie	Disk	Unm. Wheel	
SHAPE									
ISO#	065	030	125	050	120	220	160	220	
BROWNIE									
HP (12 PK)	PN	-	-	0271	0272	0273	0458	-	0453
CA (12 PK)	PN	0401	0403	0281	0282	0283	-	0405	-
CA (72 PK)	PN	H401	H403	-	-	-	-	-	-
FG (12 PK)	PN	0411	0413	-	-	-	-	-	-
FG (72 PK)	PN	H411	H413	-	-	-	-	-	-
GREENIE									
HP (12 PK)	PN	-	-	0274	0275	0276	0459	-	0454
CA (12 PK)	PN	0402	0404	0284	0285	0286	-	0406	-
CA (72 PK)	PN	H402	H404	-	-	-	-	-	-
FG (12 PK)	PN	0412	0414	-	-	-	-	-	-
FG (72 PK)	PN	H412	H414	-	-	-	-	-	-
SUPERGREENIE (yellow band)									
HP (12 PK)	PN	-	-	0277	0278	0279	-	-	454A
CA (12 PK)	PN	402A	404B	0287	0288	0289	-	-	-
FG (12 PK)	PN	412A	414B	-	-	-	-	-	-

BROWNIE®/GREENIE®/SUPERGREENIE POLISHERS

PN Product

See chart HP, CA, FG Shank (12 per pack)

KN7, PC2, WH6, Floppie, Disk

See chart CA & FG Shank (12 per pack)

Cup & Mini-Point

See chart CA & FG Shank (1/2 Gross/OneBox)

Cup & Mini-Point

See chart Unmounted Wheels (12 per pack)