



















DURA-WHITE® STONES



- Micrograined aluminum oxide grit
- Vibration-free performance
- Long lasting
- **Use Dura-White Stones for:** fine finishing on enamel, composites, glass ionomer restoratives & porcelains

SHAPE CATEGORY		CN1	CY2	FL1	FL2	FL3	IC2	IC4	KN1	KN3	RD1	RD2
SHAPE												
ISO#		024	025	030	025	030	032	045	050	070	030	024
HP (12 PK)	PN	0201	-	0202	0203	0204	-	0205	-	0206	0207	0208
CA (12 PK)	PN	0221	-	0222	0223	0224	-	0225	-	-	0227	0228
FG (12 PK)	PN	0241	0243	0240	0244	-	0245	-	0246	-	0247	0248

SHAPE CATEGORY		RE1	TC1	TC4	WH1	WH4
SHAPE						
ISO#		062	020	040	040	060
HP (12 PK)	PN	0209	0210	0211	-	0212
CA (12 PK)	PN	0229	0230	-	-	0232
FG (12 PK)	PN	-	0249	-	0250	-

DURA-WHITE® STONES
PN Product
 See chartsHP, CA, FG Shank (12 per pack)